

vive
health
portal.vivehealth.com

VIVE

AIRCARBON POWER WHEELCHAIR

MOB1140GRY

The AirCarbon Power Wheelchair is a lightweight, foldable mobility solution designed to support everyday transportation and travel needs. Constructed with a carbon fiber frame, the unit provides a balance of strength and reduced weight, allowing for easier handling, storage, and transport.

With a compact folding design and intuitive joystick control system, the wheelchair supports maneuverability in both indoor and outdoor environments. The system delivers a travel range of up to 7.5 miles per charge and operates at speeds up to 3.7 mph, making it suitable for routine mobility use. Designed to support users up to 220 lbs, the unit provides consistent performance in a portable, travel-ready format.

WHAT'S INCLUDED

- 1x AirCarbon Power Wheelchair
- 1x Battery (included).

SAFELY AND SMOOTHLY
NAVIGATING OVER ANY SURFACE



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com

DIMENSIONS:

Overall:	29" x 40" x 36" (W x L x H).
Folded Width:	24.5".
Seat:	18" x 16" x 21.25" (W x D x H).
Armrest Height:	9".
Backrest Height:	18.5".
Front Wheel Diameter:	7.75".
Rear Wheel Diameter:	9.75".
Weight:	41.9 lbs (without battery).
Weight Capacity:	220 lbs.
Maximum Speed:	3.7 mph.
Maximum Range:	Up to 7.5 miles.
Turning Radius:	47".
Maximum Incline:	6".

MATERIAL(S):

Frame:	Carbon Fiber.
Additional Materials:	Nylon.

CARE INSTRUCTIONS:

Wipe with a damp cloth as needed.
Dry thoroughly after cleaning.
Store in a clean, dry environment.

WARRANTY INFORMATION

Frame:	5 years.
Electrical:	6 months.
Batteries:	6 months.
Horn:	6 months.
Controller:	6 months.
Handles:	90 days.
Charger:	90 days.
Battery case:	60 days.
Wheels:	none / on arrival 30 days.
Screws / bolts / nuts:	none / on arrival 30 days.
Cushions / seats:	none / on arrival 30 days.
Paint / panels:	none / on arrival 30 days.
Keys:	none / on arrival 30 days.
Joysticks:	1 year.

EFFORTLESS ALL-TERRAIN MOBILITY

Navigate smoothly over carpet, tile, sidewalks, grass, and uneven paths with 9" flat-free tires and 2.5" ground clearance. Rear anti-tip wheels provide added control on inclines, making indoor and outdoor travel safer and more stable.

LIGHTWEIGHT CARBON FIBER CONSTRUCTION

Constructed from carbon fiber materials, the wheelchair maintains structural durability while reducing overall weight. This supports easier lifting, transport, and storage compared to traditional mobility devices.

FOLDABLE, TRAVEL-READY DESIGN

The frame folds quickly into a compact form factor, allowing the unit to fit in standard vehicle trunks or storage areas. This design supports portability for travel, daily errands, and limited storage environments.

COMPACT MANEUVERABILITY SYSTEM

With a 23.4" width and 47" turning radius, the wheelchair is designed for navigation through narrow spaces such as hallways, doorways, and elevators. Suitable for both indoor and outdoor use.

BATTERY-POWERED MOBILITY PERFORMANCE

The system supports a travel range of up to 7.5 miles per charge with a maximum speed of 3.7 mph. Designed for consistent daily use, the battery system provides reliable operation across typical environments.

STABLE FRAME WITH DEFINED LOAD CAPACITY

Engineered to support up to 220 lbs, the frame provides consistent stability during operation. Suitable for routine mobility applications in residential or care environments.